



SimpleReport

user guide



version 1.6 – 7/11/2022
[what's new in this version*](#)

Getting started	3
How to sign up for SimpleReport	3
SimpleReport sign up checklist	3
Fill out your organization information	3
Enter your personal contact details	3
Verify your identity	4
(optional) Get training	4
Activate your organization or user account	4
Multi-factor authentication options	5
SMS authentication	6
Google Authenticator or Okta Verify	6
Log in to SimpleReport	8
Get training	9
Using SimpleReport	9
Select your testing facility	9
Conduct and submit tests	10

Backdate a test	14
Manage results	15
Review results	15
Send someone's test results	15
Download test results	16
Correct a previous test result	17
Manage people you test	19
Self-registration	19
Add a new person	20
Archive a person	21
Update someone's profile	21
Manage users	22
Manage user permissions	22
Add new user	23
Edit user information	24
Reactivate user	25
Reset user password	26
Reset user MFA	27
Manage facility info	28
Add a testing facility or location	28
Update testing facility settings	29
Test for other diseases	31
Influenza A and B	31
Reset your SimpleReport password	31
Contact us	35
What's new in this version	35

Getting started

How to sign up for SimpleReport

If you're doing COVID-19 rapid testing, your first step is to choose just one person to sign up for your organization account. They'll create the account and complete identity verification. This person — called the account administrator — can add testing facilities and staff after creating the organization account.

What's an organization? Organizations are companies, institutions, school systems, and more, that often have multiple testing facilities or locations as part of their network. For example, John Doe School District is the "organization" that includes Highland High School, Middleton Middle School, etc. Each individual school is considered a "testing facility" that must be connected to its parent organization in SimpleReport. If you don't have a parent organization that's ok — you'll just enter in your information a bit differently when you sign up.

Let's get started.

SimpleReport sign up checklist

Complete the following steps to start using SimpleReport:

- Fill out your organization information
- Enter your personal contact details
- Verify your identity
- (optional) Get training

Fill out your organization information

You'll fill out basic information about your organization and your work contact details as part of our [SimpleReport organization sign up](#).

Make sure that the organization name that you enter in the form is accurate, since this is the name that SimpleReport will share with your public health department.

Enter your personal contact details

As the organization administrator, you'll provide your personal contact information to help verify your identity.

Verify your identity

You'll answer identity verification questions with our verification partner, Experian. This takes just a few minutes and doesn't impact your credit score.

If Experian can't verify your identity, you'll need to schedule a separate video call. During the call, you'll be asked to show two forms of ID. We recommend a driver's license and a work ID, or you can use any of the [documents on this list](#).

Your SimpleReport account will be accessible after your identity is verified. You'll simply activate and log in to your account.

(optional) Get training

Learn how to use SimpleReport:

- [Practice using SimpleReport](#) using sample data, and get familiar with different features.
- Watch these [quick training videos](#).
- Check out our step-by-step instructions for [using SimpleReport](#).


Add testing facilities and staff: Once you're all set up in SimpleReport you can [invite other staff members](#) or [add more facilities](#).

Activate your organization or user account

You'll receive an email from SimpleReport asking you to activate your account. If you're an administrator, once you've set up your account, you can [invite members of your team to SimpleReport](#).

Once you receive your account activation email, follow these steps:

1. Click the Activate your SimpleReport account button. You'll need to click the button within 7 days of receiving the invitation, or else it will expire.



2. On the page that appears, you'll need to pick:

- A secure password
- A forgot password question
- A security image

These all help to keep your account secure.

3. Click Create My Account at the bottom right of the page.

4. [Pick a multi-factor authentication option](#), and click Setup below it. These authentication options, too, are meant to secure your account. For most users, SMS authentication is the easiest option, because it sends a text message to your phone. (If you choose biometric authentication, make sure that you'll only log in to SimpleReport from the device you're currently using. Biometric authentication may be difficult to use across multiple devices.) Below, you can review instructions for setting up [SMS authentication](#) and [Google Authenticator/Okta Verify](#).

5. Click the Finish button at the bottom of the registration page. The page will confirm that you've successfully set up your multi-factor authentication, and you'll get a confirmation email, too.


You can now [log in to SimpleReport](#).

Multi-factor authentication options

If you choose SMS or Google Authenticator/Okta Verify as your multi-factor authentication, follow the instructions below to get set up.

SMS authentication


1. Enter your phone number, then click Send code. (Make sure to click Send code, or you won't be able to continue.)

A screenshot of the SimpleReport SMS authentication interface. At the top, the SimpleReport logo is displayed. Below it is a blue circular icon with a white speech bubble containing the text 'SMS'. Underneath the icon are two grey rectangular input fields. Below these is a 'Phone number' label followed by a small grey box containing '+1' and a larger white input field. To the right of the input field is a blue button labeled 'Send code'. At the bottom of the form is another grey rectangular input field.

2. Check your text messages for a 6-digit authentication code.
3. Enter the code in the Enter Code field, then click Verify.

Google Authenticator or Okta Verify

1. Select the kind of phone that you use (either iPhone or Android). You'll be asked to download an app. Download it on your phone and wait for it to install. (The page on your device might vary a bit from the screenshot below, based on your device type and whether you chose Google Authenticator or Okta Verify.)



2. Once the app is installed, click **Next**.
3. Open the app and scan the QR code that appears on your SimpleReport registration page. Once you've successfully scanned the QR code, click **Next**.



4. Back on your phone, the app will show you a code. Enter the code on the registration page, then click **Verify**. (The code changes regularly, so you'll need to check the app each time you log in to SimpleReport.)

Log in to SimpleReport

Once you've activated your account, you can log in and start using SimpleReport.

To log in:

1. Go to SimpleReport.gov and click Log in at the top right corner of the page.
2. Enter your username and password.
3. You'll need to use the multi-factor authentication option that you set up when you registered your account. Here's an example of what you'll need to do if you set up SMS authentication:
 - a. Click the gray **Send code** button to have the verification code sent to your phone. (Make sure to click the button, or else the code won't get sent.)



- b. Once you receive the text, enter the code in the Enter Code field.
 - c. Click **Verify**.
4. If your organization has multiple testing facilities, select yours from the choices shown.

Once you've successfully logged in, you can do a number of tasks, like [conduct a test](#) or [add a new person to SimpleReport](#).

Get training


Check out our resources to get familiar with SimpleReport's testing, reporting, and workflow features.

- Check out our video [introduction and onboarding guide](#).
- Practice using SimpleReport and take a look around on our [demo site](#).
- Have a look at our [resources page](#).
- Review our [K-12 schools guide](#).
- Still have a question? Visit our [support page](#).

Using SimpleReport

Select your testing facility


If you work at multiple testing facilities, you'll be asked about the testing facility where you're working each time you log in to SimpleReport. Just choose the correct testing facility and you'll be ready to [conduct tests](#).




The image shows a screenshot of the SimpleReport login interface. At the top, the 'SimpleReport' logo is displayed in a dark blue font. Below the logo is a grey rectangular placeholder. Further down is another grey rectangular placeholder. Below that, the text 'Please select which facility you are working at today' is centered. Underneath this text are three rectangular buttons with blue borders and blue text, stacked vertically. The buttons are labeled 'Your current facility', 'Another facility', and 'One more facility'.

You might have already chosen a testing facility, and need to change it. To switch between testing facilities after you've already chosen one:

1. At the top right corner of the page, just to the left of the gear and person icons, find the name of your current testing facility. Click it to reveal a dropdown menu of other facilities.




2. Choose the facility you want to switch to.




Conduct and submit tests

Once you've [logged in to SimpleReport](#), you can begin conducting tests and submitting results to your public health department:

1. Make sure that you're on the "Conduct tests" page. It's the page that appears once you've logged in and [picked the facility where you're working](#). You can also get to it by clicking **Conduct tests** at the top of the page.




2. In the text field just below **Conduct tests**, enter the name of the person you're testing.



The image shows the 'SimpleReport' logo on the left and a dark blue header with the text 'Conduct Tests' on the right. Below the header is a light gray search bar with the placeholder text 'Search for a person to start their test' and a magnifying glass icon on the right.

3. Confirm the person's full name and date of birth. Then click **Begin test** to the right of their name.




The image shows a search results dropdown menu. At the top, there is a search bar with the text 'jane' and a blue search button with a magnifying glass icon. Below the search bar is a table with three columns: 'Full name', 'Date of birth', and 'Actions'.

Full name	Date of birth	Actions
Jane Doe	1980-01-10	Begin test
Jane Abbey Doe	1996-06-28	Begin test

If the person you're looking for doesn't appear in the search results, make sure that you've spelled their name correctly. If you still can't find them, you'll need to [add them](#) before proceeding with the rest of these steps.

4. The person you're testing will need to answer a test questionnaire.



The image shows a digital form interface. At the top, there is a section labeled 'Results' which contains several horizontal bars of varying lengths, likely representing a scale or checklist. Below this is another section labeled 'Symptoms' which also contains horizontal bars. At the bottom right of the form, there is a blue button labeled 'Continue'.

The first question asks the patient how they want to get their result. If they choose to get their result by text or email, they'll get a notification as soon as you submit their result in SimpleReport. (The text and email messages include a link to a secure portal to check their result.) If they don't want their results via text or email, you can also [print a copy](#). The questionnaire also has a section to enter information about any symptoms the patient may have at the time of the test.

You can ask these test questions now, or click **Continue** and come back to them later. Either way, you'll need to complete them before submitting test results.

5. The person's name will appear on a card on the "Conduct tests" page. You can continue with their test or, if you're testing multiple people at once, you can repeat the previous steps for each of the people you're testing.
6. Collect the swab from the person you're testing and prepare their sample.
7. *(optional)* On the person's card, click the blue button with the stopwatch icon to start the timer. If your device's sound is on, you'll hear an alarm at the end of the required time. You can restart the timer at

any time by clicking the blue button again.

A screenshot of a user profile card for Jane Doe. The card has a header with the name "Jane Doe" on the left and a blue button with "15:00" and a clock icon on the right. Below the header, there are three rows of placeholder text represented by gray rectangles.

8. You'll have the option to select the testing device and swab type from the dropdowns. Make sure you select the correct device and swab type. (It defaults to the last device used.) If the device or swab type you're using isn't on the list, ask an administrator to [add the device or swab type to your SimpleReport facility settings](#). If the device isn't available for them to add, the administrator will need to contact SimpleReport to request this.

A screenshot of the testing device and swab type selection interface. It shows a form with a "Device" dropdown menu and a "Swab type" dropdown menu. The "Device" dropdown is highlighted with a blue border. To the right of these dropdowns, there is a text input field with a placeholder gray rectangle.

9. Once the test is complete, select the result of the COVID-19 test on the right side of the person's card: "Positive", "Negative", or "Inconclusive".
10. If the device you're using also tests for flu (a "multiplex device"), you'll see columns for flu A and B results next to the COVID-19 column. You must enter "Positive" or "Negative" for each illness, or select "Inconclusive" for all three before you can submit.

A screenshot of the test result selection interface. It shows a form with a header for "Doe, Jane" and a close button (X) in the top right corner. Below the header, there are three columns for "COVID-19", "Flu A", and "Flu B". Each column has two radio buttons: "Positive (+)" and "Negative (-)". The "Negative (-)" radio button for COVID-19 is selected. Below the radio buttons, there is a checkbox labeled "Mark test as inconclusive". To the right of the checkbox, there is a blue button labeled "Submit" with an information icon (i) to its left.

11. Make sure that you've completed the test questionnaire. You can review the questions by clicking **Test questionnaire** on the card. A gray "PENDING" tag means that some answers are still missing. A green "COMPLETE" tag means that all questions are completed.

You can submit the results without completing all of the questions, if the person being tested declines to answer them.


12. Click **Submit**.

You've successfully submitted the test result to your public health department. You're all done. If you have the SimpleReport permissions to do so, you can also [review submitted results](#).

Backdate a test

If you're submitting test results you've already conducted, you can backdate the result in SimpleReport.

1. Go to **Conduct tests**, search for a patient name, and click **Begin test**
2. Uncheck the checkbox that says **Use current date**

A screenshot of the SimpleReport web application interface. At the top, there is a navigation bar with the SimpleReport logo and a 'Conduct Tests' button. Below the navigation bar is a search bar with the placeholder text 'Search for a person to start their test'. The main content area is divided into two columns. The left column contains a form for entering patient information, including fields for name, date of birth, and gender. The right column contains a form for entering test results, including a dropdown for 'Test type' and a text area for 'Test result'. In the bottom right corner of the form, there is a checkbox labeled 'Use current date' which is currently checked. A mouse cursor is pointing at the checkbox, indicating it should be unchecked.

3. Enter the correct Test date and Test time. Enter the result by clicking **Submit**.

Test date
mm/dd/yyyy

Test time
hh:mm


Submit

Manage results

Review results

If you work for a testing facility, here's how to find your patient's COVID-19 results:

1. Click **Results** at the top of the page.



2. A table of results from your testing facility will appear showing all results from the past 48 hours. To show older results, you can always click the blue See all results button at the bottom of the page.


If there's a problem with any result, you can [mark it as an error](#).

Send someone's test results

You can send or resend test results to patients through SimpleReport.

To send or resend someone's test result:


1. Click **Results** at the top of the page.



2. A table of results from your testing facility will appear showing all results from the past 48 hours. To show older results, you can always click the blue **See all results** button at the bottom of the page.
3. Find the result that you want to send, and go all the way to the right of the page. Under the “Actions” column, click the three dots in the same row as the result that you want to send.

Test Results (all)				
				Actions
Jane Doe				...

4. Select the option you want then select **Print result**, **Email result**, or **Text result**.



5. A preview of the test result will appear. Select the option you want and follow the instructions to print, email, or text the result.

Download test results

You can filter results by name, test date, test result, or role and download them as a CSV file to your computer.

1. To download test results:
2. Enter the information to filter and select the type of test results you want to download. Filters are any of the categories at the top of the page, for example a test date or range, a type of result, and/or results from a certain test group (role) such as staff, students, visitors, etc.

Test results
Showing 1-2 of 2

Download results
Clear filters

Search by name

Date range (start)

Date range (end)

Result

Role

-Select-

-Select-

3. Select **Download results** and confirm. The results with the filters you selected will download to your computer.
4. To view other test results or start a new download, select **Clear filters**.

Correct a previous test result

If something goes wrong with a test result, you can correct it in SimpleReport. The previous result will still be visible, but marked as an error and submitted to the public health department. Make sure to also submit a new, corrected test result (backdated to the original test date) after completing the steps here.

To correct a previous test result:

1. Click **Results** at the top of the page.



2. A table of results from your testing facility will appear showing all results from the past 48 hours. To show older results, you can always click the blue **See all results** button at the bottom of the page.
3. Find the result in the table, and go all the way to the right of the page. Under the "Actions" column, click the three dots in the same row as the result that you want to correct.

Test Results (all)				
				Actions
Jane Doe				...

- Click **Mark as error**.



- A box will appear and ask you why you want to mark the result as an error. Enter a reason, then click **Yes, I'm sure**. (Your reason needs to be at least 4 characters long before you can proceed.)

A confirmation dialog box with a dark border. It contains the text: 'Are you sure you want to mark this test result for **Jane Doe** as an error?'. Below this is the text: 'If so, please enter a reason.' followed by a text input field. Inside the input field, there is a placeholder text: 'Make sure that your reason is at least 4 characters long.' At the bottom of the dialog, there are two buttons: a blue link-style button labeled 'No, go back' and a solid blue button labeled 'Yes, I'm sure'.

The result will be marked as an error, and will appear with a strike through it when you review it in the future.

Test Results (all)				
Person Name	Date of Test	Result	Device	Actions
Jane Doe	Feb 12, 2021 11:54 AM	NEGATIVE	Abbott IDNow	

To submit a new, corrected result to replace the one now marked as an error, follow the instructions to [conduct and submit a test](#), and be sure to update the date and time of the test to reflect when the test originally occurred (you may need to backdate it).

Manage people you test

Self-registration

The fastest way to get new people into SimpleReport and start conducting tests is through patient self-registration. Only new patients need to self-register.

To offer self-registration to your patients:

1. Go to the gear icon, then click on the “Patient self-registration page” on SimpleReport. Use the “Organization link” if you want the patient to be able to get tested at any of your locations, or use the relevant “Testing facility link” if they will only be tested at one location.

Testing facility name	Patient self-registration link
Testing Site	training.simplereport.gov/register/T37FT

2. Send the link to new patients before their test via email, text message, mail, etc. You can also give out the link in person at your testing facility using a sign or printed form with a QR code, for example.
3. Patients complete the self-registration form online, at their convenience.
4. When they arrive for testing, simply search for patients by name in SimpleReport and follow your normal process to check them in.


Add a new person

In order to report test results for someone, you'll need to make sure that their information has been added to SimpleReport.


You can find out if they've already been added by searching for their name in the "Search for a person to start their test" field, just beneath the SimpleReport logo on the "Conduct Tests" page. If their name doesn't appear in the search results, you can add them.

To add someone new:

1. Click **People** at the top of the page.



2. Click the **+ Add person** button on the right side of the page.



3. Enter the person's information on the "Add New Person" page that appears. There are a few types of information to enter for the new person:
 - **Role (optional):** Mark the person as either "Staff", "Resident", "Visitor", or "Student".
 - **Facility (required):** Associate the person with all of your organization's testing facilities, or just one. They'll only appear in search results for tests conducted at the facilities that they're associated with.
 - **Demographics (optional, but encouraged):** Ask the person how they identify in terms of their race, ethnicity, and sex. We know that public health problems are disproportionately high in some populations. Filling out this information can help efforts to recognize and mitigate disparities in health outcomes.

Many questions here appear as options. Please refer to your public health department requirements to make sure you complete all required information.

4. When you're done, click **Save changes** at either the bottom right or top right corner of the page.

A green confirmation box will appear to let you know that you've successfully added someone new. You're all done.

If you're ready to conduct a test for this person, learn how to [conduct and submit tests](#).

Archive a person

Archiving a person removes them from your list of people. Because it's archived and not deleted, their record will still be accessible if it needs to be recovered later. (To unarchive a person, you'll need to [contact SimpleReport support](#).)

To archive a person:

1. Click **People** at the top of the page.




2. Go to the Actions column.
3. Click on the three dots next to the person's name you want to archive, then click **Archive record**.
4. Confirm archive request.

Update someone's profile

After you've added someone in SimpleReport, you may need to return and update their information, including adding or modifying the testing facilities where they'll be tested.

To update someone's profile:

1. Click **People** at the top of the page.



2. Click on the name of the person whose profile you want to update.
3. A page with all of the person's information will appear. Add or edit whichever information you need to update.

4. Click the blue **Save changes** button at either the top right or bottom left of the page. (The button will be gray until you've made a change.)

You'll be sent back to the "People" page, and a green box will appear at the bottom of the page to confirm that SimpleReport has saved your update. You're all done.

Manage users


Manage user permissions

You might want to limit which members of your staff can perform certain tasks or access certain information in SimpleReport. All users can conduct and report tests, but you can restrict access to other things, like viewing and correcting results, viewing and editing profiles, and accessing facility settings.


To manage user permissions in SimpleReport, you must have admin privileges. If you don't have admin privileges, talk to an admin about upgrading permissions on your account.

To change permissions for a user:

1. Click the gear icon at the top right corner of the page.



2. Beneath the SimpleReport logo at the top of the page, find the "Manage users" tab. Click "Facility access" below.




3. Under "User role", click the user whose access you want to edit. There are three levels of access:

“Admin”, “Standard”, and “Testing only”:

- **Admin** users have full access to conduct and report tests, manage test results, manage patient profiles, and manage testing facility and user (staff) settings.
- **Standard** users can conduct and report tests, manage test results, and manage patient profiles. They can’t manage testing facility or user (staff) settings.
- **Testing only** users can only conduct and report tests.

4. Select the access level that you want to give this user, then click **Save changes**.



Manage users

Facility access

User role

Admins have full access to use and change settings on SimpleReport. Standard and testing-only users have limited access for specific tasks, as described below.

☐ **Admin**
Full access: Conduct tests, manage test results and patient profiles, manage account settings, users, and testing facilities.

☒ **Standard user**
Conduct tests, manage test results, and patient profiles


☐ **Testing only**
Conduct tests

The next time this user logs in to SimpleReport, they’ll have the level of access that you just set.


Add new user

Admins can give staff and employees SimpleReport access by adding them directly in SimpleReport. To add a new user:

1. Click the gear icon at the top right corner of the page.



2. Beneath the SimpleReport logo at the top of the page, find the “Manage Organization”, “Manage Facilities”, and “Manage Users” tabs. Click **Manage Users**.




3. Under “Users”, click **+ New user** in the top right corner.
4. Enter the first name, last name, and email address of the person you want to add, then click **Send invite**. They’ll receive an email with a link to sign up for a SimpleReport account.
5. New users will automatically be assigned the “Standard” access level, which allows them to conduct tests and manage test results and patient profiles. Once you send the invite, you can change the user’s access level by following the instructions to [manage user permissions](#).

Edit user information


Admins can update or edit basic user information such as name and email for staff and employees.

To edit user information:

1. Click the gear icon at the top right corner of the page.



2. Find and click on the “**Manage users**” tab.



3. Select **Edit name** or **Edit email** and confirm your changes.

Reactivate user


SimpleReport deactivates users after 60 days without logging in. Organization admins can reactivate users in SimpleReport.

To reactivate a user:

1. Click the gear icon at the top right corner of the page.



2. Find and click **Manage users**.



3. Find the user that needs to be reactivated and click the **Activate user** button on their profile and confirm.


Please note: If the reactivated user doesn't log in to SimpleReport before 6AM Eastern time the following day, their account will be deactivated again.

Reset user password


Organization admins can reset passwords for staff or employees (users) in SimpleReport. If you don't have admin privileges, talk to an admin about upgrading permissions on your account. (You can also [reset your own password](#).)

To reset a user password:


1. Click the gear icon at the top right corner of the page.



2. Find and click "Manage users".



3. Select **Send password reset email**, and confirm.




Reset user MFA


Organization admins can reset multi-factor authentication (MFA) login access settings for staff or employees (users) in SimpleReport. If you don't have admin privileges, talk to an admin about upgrading permissions on your account.

To reset a user MFA:


1. Click the gear icon at the top right corner of the page.



2. Find and click "Manage users".



3. Find and click on the correct user from the list, select **Reset MFA**, then confirm.




Manage facility info

Add a testing facility or location


You can add a testing facility or location directly in SimpleReport. Before you add a new facility, make sure that it's in [a jurisdiction that SimpleReport supports](#).

To add a testing facility:


1. Click the gear icon at the top right corner of the page. (If the gear doesn't appear in your account, you'll need to have an administrator [change your account permissions](#) before you can add a facility.)




2. Beneath the SimpleReport logo at the top of the page, find the "Manage organization" and "Manage facilities" tabs. (Your organization is the umbrella for all of your facilities.) Click **Manage facilities**.



3. Click **+ New facility**.



4. Enter information for the new facility. This includes facility contact information, CLIA number, ordering provider contact information and NPI number, and test devices. (Links to look up CLIA and NPI numbers are provided on this page.)
5. When you're done, scroll back to the top right of the page and click **Save changes**.




Update testing facility settings

You may need to change information about your testing facility, from the phone number to the ordering provider to the testing devices you use.


Only change testing facility names if doing so is absolutely necessary. To change the name of your facility, contact support@simplereport.gov to make sure that your results continue being reported to your public health department.

To update your testing facility information:


1. Click the gear icon at the top right corner of the page.



2. Beneath the SimpleReport logo at the top of the page, find the “Manage Organization” and “Manage Facilities” tabs. (Your organization is the umbrella for all of your testing facilities.) Click **Manage Facilities**.



3. Click the name of the testing facility whose information you want to edit. (You can also [add a new facility](#).)



4. Edit the information that you want to change.

If the testing devices your facility uses have changed, update your device information under “Manage devices” card.



Device type	Action
Cue	
Cepheid Xpert Xpress	

[+ Add device](#)

5. When you're done, scroll back to the top right of the page and click **Save changes**.

Test for other diseases

Influenza A and B

Certain devices test for COVID-19, influenza A, and influenza B (“multiplex devices”). SimpleReport requires you to include flu results when you record tests from these devices. Flu results are not sent to public health departments, but you can share them with patients.

To record flu results on SimpleReport, [conduct and submit a test](#) as you normally would. You'll be able to enter results for flu A and flu B, along with COVID-19, when you select a multiplex device.


Devices that test for flu and COVID-19 include:

- CDC Influenza SARS-CoV-2 Multiplex Assay (RT-PCR)
- Abbott Alinity M (RT-PCR)
- BD Veritor System for Rapid Detection of SARS-CoV-2 & Flu A+B (Antigen)
- Quidel Sofia 2 (Antigen)

Reset your SimpleReport password


To reset your password you can follow these steps:

1. Go to SimpleReport.gov
2. Click **Log in** at the top right of the page.
3. Click **Need help signing in?** at the bottom of the page.
4. Click **Forgot password.**
5. Enter your email address, then click the blue **Reset via email** button.




The screenshot shows the 'Reset Password' page on SimpleReport.gov. At the top is the SimpleReport logo. Below it is the title 'Reset Password'. There is a text input field labeled 'Email or Username'. Below the input field are two horizontal gray bars. Then there are three blue buttons: the first has a white horizontal line, the second also has a white horizontal line, and the third is labeled 'Reset via Email'. At the bottom is another horizontal gray bar.

6. Check your email inbox. You should receive an email allowing you to reset your password. Open the email and click the link inside.
7. Answer your security question, then click **Reset Password.**




The screenshot shows the 'Answer Forgotten Password Challenge' page on SimpleReport.gov. At the top is the SimpleReport logo. Below it is a circular placeholder for a profile picture. The title 'Answer Forgotten Password Challenge' is centered. Below the title is a horizontal gray bar, followed by a text input field labeled 'Answer'. Below the input field is another horizontal gray bar. Then there is a blue button labeled 'Reset Password'. At the bottom is a final horizontal gray bar.

8. Enter a new password in both the “New password” and “Repeat password” fields, then click **Reset Password**.




The screenshot shows a web form titled "Reset your Okta password". At the top, there is a list of bulleted items, each followed by a gray rectangular placeholder. Below this list are two input fields: "New password" and "Repeat password". At the bottom of the form is a blue button labeled "Reset Password".

9. You'll need to use the multi-factor authentication option that you set up when you registered your account. Here's an example of what you'll need to do if you set up SMS authentication:
 - a. Click the gray **Send code** button to have the verification code sent to your phone. (Make sure to click the button, or else the code won't get sent.)




The screenshot shows the SimpleReport login interface. At the top is the SimpleReport logo. Below it is a circular icon with a blue speech bubble and the text "SMS". Underneath are several gray rectangular placeholders. There is an input field followed by a gray button labeled "Send code". Below that are more gray placeholders and a large blue button with a white checkmark icon.

- b. Once you receive the text, enter the code in the “Enter Code” field.

A screenshot of the SimpleReport SMS verification interface. At the top is the SimpleReport logo. Below it is a blue circular icon with a white speech bubble and the text "SMS". There are two horizontal bars for text input. Below these is the label "Enter Code" followed by a text input field and a "Sent" button. At the bottom is a large blue button with a white checkmark icon.

- c. Click **Verify**.

10. Once you’ve finished verifying, the Okta page will appear. Click **Simple Report**.



11. You’ll be sent back to the SimpleReport landing page. Click **Log in** at the top right corner of the page, which should automatically log you back into SimpleReport.

If you’re still having trouble changing your password, SimpleReport admins can also [reset user passwords](#) for testing facility employees and staff.

Contact us

If you've already read our [troubleshooting techniques](#) and are still having problems with SimpleReport, please reach out to us.

Contact support@simplereport.gov with information about the problems you're experiencing. Tell us what's going on and what you're trying to do, and include screenshots of the issue. We'll follow up to help.

What's new in this version

Version 1.6 (7/11/22) of the SimpleReport user guide includes:

- An additional step in "[Conduct and submit tests](#)" for devices that also test for flu
- A new section about testing for other diseases, with information about [Influenza A and B](#)
- A new section with instructions to [reset your own SimpleReport password](#)
- A link in the "Reset user password" section to guidance on resetting your own password